

How to use the ringsize meter:

- 1 - Read the contents of the page carefully before printing.
- 2 - Cut out the ringsize meter and make sure to also cut the slit on the large part of the meter (-----O).
- 3 - Slide the pointed tip into the slit making sure that the coloured side is facing outwards.
- 4 - Slide the ringsize meter over the desired finger at the first phalanx and pull the pointed tip until it fits around.
- 5 - Read the number visible at the slit - which corresponds to your finger's circumference in millimeters - and continue to page 2, if you wish.
- 6 - The table on page 2 features circles with two numbers:
 - a) the number inside the circle corresponds to the circumference
 - b) the number underneath the circle indicates its diameter.You can use this page to identify the size of a simple ring band laying it flat over the corresponding circle.

Note: The achieved measures are indicative. KB can not be held responsible for possible measuring mistakes.



In order to verify that the ringsize meter was printed correctly, please measure with a ruler the design to the right. It should be exactly 2 centimeters long.





40

12,69mm



41

13,05mm



42

13,37mm



43

13,69mm



44

14,09mm



45

14,33mm



46

14,64mm



47

14,96mm



48

15,28mm



49

15,6mm



50

15,92mm



51

16,24mm



52

16,56mm



53

16,87mm



54

17,19mm



55

17,51mm



56

17,83mm



57

18,15mm



58

18,47mm



59

18,87mm



60

19,1mm



61

19,42mm




62

19,74mm



63

20,06mm



64

20,38mm



65

20,7mm



66

21,02mm



67

21,33mm



68

21,65mm



69

21,97mm



70

22,29mm